Technical Information Memobase Plus CYZ71D

Multichannel multiparameter software for liquid analysis with digital Memosens sensors



Measure, calibrate and document your Memosens sensors with one single tool

Application

- The Memobase Plus manages the complete life cycle of pH, ORP, conductivity and oxygen sensors with robust Memosens technology.
- The software can be used in all industries and meets the highest demands of the pharmaceutical industry.

Your benefits

- Greater efficiency with easy sensor maintenance
- Advanced diagnostics with "Ready for next batch" indication
- Better process safety thanks to sensor traceability
- Full flexibility with multichannel and multiparameter functionality
- 100 % consistency between lab and process measurements
- Highest accuracy for your measurement values
- Easy buffer management

Detailed information on the product benefits is available on the product page: www.endress.com/cyz71d



Function and system design


Measuring system

A complete measuring system consists of:

- PC or Windows-based tablet with Memobase Plus software package installed and database connected
- MemoLink sensor terminal box (connection to PC, Ex barrier)
- Thin, flexible CYK20 Memosens laboratory cable or CYK10 Memosens process cable
- USB cable to connect the MemoLink sensor terminal box and PC
- Memosens sensor
- i

A PC or Windows-based tablet is not included in the delivery.

Memosens sensors must be ordered separately. Information on this can be found at: www.endress.com/memosens



A00316

- 1 Measuring system for Memobase Plus CYZ71D
- 1 PC (not supplied)
- 2 USB hub (optional, not supplied)
- 3 1 to 4 USB cables
- 4 1 to 4 MemoLink sensor terminal boxes
- 1 to 4 CYK20 Memosens laboratory cables or CYK10 Memosens process cables
- 6 1 to 4 Memosens sensors

Connection

- USB → MemoLink sensor terminal box to PC
- Memosens data cable → Sensor to MemoLink sensor terminal box

System requirements

System requirements for installation and use of Memobase Plus:

System requirements

Operating system Windows 7 Service Pack 1 (32 and 64 bit)

Windows 10 (32 and 64 bit)

Monitor At least 1280×1024 pixel, also suitable for touchscreen

Processor Minimum clock speed of 1 GHz

Not a virtual machine

RAM 1

USB At least a type A USB interface

At least USB 2.0

Barcode reader Supported interfaces:

USB-HID interfaceUSB-COM interface



The interface must be configured on the barcode reader.

Minimum resolution: 0.254 mm (10.0 mil)

Other $\hfill \mbox{ \ \, } \mbox{ \ \, } \mbox{ CD-/DVD drive or internet access for program installation}$

Adobe ReaderPrinter driver


Microsoft .NET Framework 4.7

Software functions

Memobase Plus has four main functions, which are listed in the navigation bar on the left:

- Measure: measurement including graph and description of sample
- Calibrate: several calibration methods and testing equipment management
- Sensors: settings, administration, status and information
- Reports: database view, report creation and export function

A separate tab is displayed at the top for each sensor connected to a MemoLink sensor terminal box. The tab displays the sensor type, order root, serial number and tag name.



■ 2 Program structure

- 1 Tabs
- 2 Main menu
- 3 Instruction area
- 4 Status bar

Measurement

- Numerical and graphical display of primary and secondary measured values (with zoom function and time bar)
- Sample description for the verifiable assignment of a measurement
- Information on measurement settings is displayed in order to avoid interpretation errors

"Advanced diagnostics" license:

- Sample data easily transmitted by scanning the barcode ¹⁾
- Measuring range monitoring

The sensor measuring range is highlighted in color in the graph:

- Measured values within the sensor measuring range meet the GLP requirements and can be exported and saved
- Measured values outside the sensor measuring range cannot be exported and saved




Measuring range monitoring ("Advanced diagnostics" license): the measuring range including tolerances is highlighted in color

¹⁾ Prerequisite: generated barcode contains the relevant data (for detailed information on barcode specifications, see Operating Instructions BA00502C)

Calibration and adjustment

- Guided step-by-step calibration with clear instructions
- Reference solution management with preprogrammed values for the most common buffer solutions (pH) available on the market
- Live-graph for visual monitoring during calibration enables sensor condition appraisal
- Ability to adapt stability criteria to different requirements for optimized measuring performance
- Optional "as-found-as-left" report provides important information regarding the sensor performance and the consistency of the current process



Multipoint calibration ("Advanced diagnostics" license): buffers 4 and 7 selected as the adjustment points, buffer 9 as the calibration point

"Advanced diagnostics" license

- Multipoint calibration and adjustment with up to 10 measured values from pH sensors
- For pH: monitoring limits can be defined for deviation between measured value and known standard
- Easy transfer of data from Endress+Hauser testing equipment by simply scanning the barcode ²⁾
 - pH: CPY20 buffer solutions ²⁾
 - Conductivity: CLY11 calibration solutions ²⁾
 - Oxygen: COY8 zero point gel²⁾
- Memobase Plus supports user administration functions, electronic documentation and signatures in accordance with "Food and Drug Administration (FDA)" 21 CFR Part 11

The complete range of functions of the audit trail is available only with the "Pharmaceutical compliance" license.

The "Memobase Plus Basic" and "Advanced Diagnostics" licenses allow read-only access to the diagnostic messages in the audit trail.


Sensors

- Time stamp for adjustment and deactivation (with explanation) documents the entire sensor life cycle
- Possible to assign sensor to measuring point in the lab
- Specification of calibration method for effective, schedulable work
- Operating hours counter to analyze the sensor condition
- Validity of sensor calibration checked on an hourly basis ("Advanced diagnostics" license)
 - Define intervals on an hourly basis for the calibration and adjustment of sensors
 - Alarms and warnings alert users to pending calibration and adjustments

²⁾ Prerequisite: current Endress+Hauser testing equipment with relevant barcode

Reports

- Subdivision into Measure / Calibrate / Sensors / Test equipment categories and categorization by measuring parameter enables the fast retrieval of data
- Sorting and filter function helps users find data more quickly in every column
- Reports at the touch of a button, optionally with company's own logo
- The report contains all of the required information, including a table with new calibration values, deviations from old values as well as calibration history charts (slope and zero point)
- Export to .PDF, .XML- or .CSV file for further processing and analysis, e.g. in Microsoft Excel or LIMS systems



■ 5 Report creation

setur

- Audit trail compliant with the requirements of the pharmaceutical industry, and user administration with five roles for full traceability
- Languages:
 - German
 - English
 - Spanish
 - Italian
 - French
 - Dutch
 - Portuguese
 - Polish
 - Czech
 - Russian
 - Turkish
 - Japanese
 - Chinese
- Database settings including test function and initialization

Diagnostic messages

- Diagnostic messages are characterized as per Namur NE 107, including the corresponding symbols
- A window with instructions on how to proceed appears when quality- and safety-related messages are displayed
- All other messages are displayed in the status bar

Network architecture

Memobase Plus is based on a client-server architecture and allows several clients to access a shared central database.

Supported databases:


- Microsoft SQL Server (included in delivery)
- Oracle (interface available)


Possible installations:

- Local installation on a PC or Windows-based tablet
- Central installation for simultaneous use by multiple PCs or Windows-based tablets

Possible operating modes:

- Master mode:
 - A local or central database is connected to Memobase Plus
- Master-slave mode:
 - A central database is set up as the "master" and one or more local databases act as the "slave"
 - Data can be saved in a local database and transmitted to a central database at a later stage






A0031877

₽ 7

■ 6 Example 1 for master mode: installation with 1 central database which 4 clients access

Example 2 for master mode: installation with 1 local database



A003187

- \blacksquare 8 Example for master-slave mode: installation with 1 central and 4 local databases
- 1 Central database (master)
- 2 Local databases (slaves)

Operation with local and central database (master-slave mode)

Enjoy complete mobility with Memobase Plus:

- Save measurement and calibration data on your PC or Windows-based tablet to a local database.
- The next time you connect to the network, the values and sensor data saved locally can be easily replicated with a central database.

Replicated data:

Master to slave

- Templates
- Testing equipment
- Specifications from user administration

Slave to master

- Sensor data
- Measuring and calibration data
- Testing equipment recorded in slave database
- Data recorded in the audit trail ("Pharmaceutics compliance" license)

Memosens technology

Memosens MEMO(SENS

Memosens makes your measuring point safer and more reliable:

- Non-contact, digital signal transmission enables optimum galvanic isolation
- No contact corrosion
- Completely watertight
 - Can even be connected under water
 - No contact corrosion
- Sensor can be calibrated in a lab, thus increasing the availability of the measuring point in the process
- Predictive maintenance thanks to recording of sensor data, e.g.:
 - Total hours of operation
 - Hours of operation with very high or very low measured values
 - Hours of operation at high temperatures
 - Number of steam sterilizations
 - Sensor condition

MemoLink input

Type of input

Memosens port: M12 socket

Measured values

All sensors with an inductive Memosens plug-in head are suitable for connection (pH/ORP, conductive conductivity and dissolved oxygen) and inductive conductivity with a fixed cable and M12 connector. All sensors contain a temperature sensor.



For detailed information on "Measured variables", see the Operating Instructions for the connected sensor.

MemoLink output

Output type

- USB port: mini USB 2.0 Type B
- USB class: HID

Output voltage

2.8 to 3.3 V

Output current


10 mA

Power supply

Supply voltage

The PC powers the sensor(s) and the MemoLink sensor terminal box(es) via the USB cable and enables the bidirectional transfer of Memosens data. If a USB hub is used, it must have a power unit.

Connection



A003165

- 1 Cable with mini USB plug
- 2 Cable with M12 plug

Power connection

- 5 V DC via USB
- Low power mode: max. 100 mA as per USB specification 2.0

Cable length

- USB cable: 2.0 m (6.6 ft)
- Memosens lab cable CYK20: ¹.5/3.0 m (⁴.9/9.8 ft) (depending on order version)
- Memosens process cable CYK10: 3 to 100 m (9.8 to 328.1 ft) (depending on order version)

Performance characteristics

No corrupted measured data



For detailed information on "Measured error", see the documentation for the connected sensor.

MemoLink only transmits data digitally so no measured data can be corrupted. The measuring signal is converted to digital data in the sensor, which means that the measured values are not affected by MemoLink, the cable or the software.

Environment

Ambient temperature range

- MemoLink: -10 to 50 °C (14 to 122 °F)
- Memosens lab cable CYK20: -10 to 50 °C (14 to 122 °F)
- Memosens process cable CYK10: -25 to 135 °C (-13 to 277 °F)

Storage temperature

- MemoLink: -25 to 85 °C (-13 to 185 °F)
- Memosens lab cable CYK20: -10 to 50 °C (14 to 122 °F)
- Memosens process cable CYK10: -25 to 135 °C (-13 to 277 °F)

Relative humidity

Maximum 85 %, non-condensating

Degree of protection


- MemoLink: IP 65 (mated, i.e. when cables are connected) in accordance with EN 60529 and Type 2 in accordance with UL
- CYK20 Memosens laboratory cable: IP 68
- CYK10 Memosens process cable: IP 68

Electromagnetic compatibility

Interference emission and interference immunity as per EN 61326-1:2006, Class B (Industrial)

Mechanical construction

Dimensions



A003165

■ 9 Dimensions of MemoLink in mm (in)



The MemoLink sensor terminal boxes can be stacked on top of one another. In such situations, the "Power / Data" LED is still easily visible.

Weight

0.24 kg (0.53 lb.) not including cable

Materials

- Housing: PBT
- Housing feet: EPDM

Certificates and approvals

C€ mark

The product meets the requirements of the harmonized European standards. As such, it complies with the legal specifications of the EU directives. The manufacturer confirms successful testing of the product by affixing to it the CC mark.

Ex approval

- MemoLink: ATEX II (2) G [Ex ia Gb] II C
- EMC Directive 2004/108/EC



The measuring point may be operated only in non-hazardous areas. Memosens sensors with and without Ex approval may be connected alternately to the Memosens interface. Connecting Memosens sensors without Ex approval does not affect the intrinsic safety of any Ex-rated Memosens sensors connected at a later stage.

Background: ATEX-certified instruments formally lose their approval as soon as they are connected to non-certified equipment. MemoLink has been developed and certified in a way that prevents this.

Ordering information

Ordering informatio

Product Configurator

Product page

On the product page there is a **Configure** button to the right of the product image.

1. Click this button.

www.endress.com/cyz71d


- └ The Configurator opens in a separate window.
- 2. Select all the options to configure the device in line with your requirements.
 - In this way, you receive a valid and complete order code for the device.
- 3. Export the order code as a PDF or Excel file. To do so, click the appropriate button on the right above the selection window.
- For many products you also have the option of downloading CAD or 2D drawings of the selected product version. Click the **CAD** tab for this and select the desired file type using picklists.

Licensing model

One license is required per workstation. The license can be used for an unlimited period of time. It is connected to the PC or Windows-based tablet that was used to generate the activation code for connection

 $1\ {
m to}\ 4\ {
m MemoLink}\ {
m sensor}\ {
m terminal}\ {
m boxes},$ with the same number of Memosens sensors (1 to 4), can be connected per license.

The licenses can be ordered as single licenses or as a multi-user license with 2 to 5 licenses.



License function scope

The functional range depends on the order configuration.

The following functional packages are available:

License	Function range
Memobase Plus Basic	Measure, calibrate, document
Advanced diagnostics	Functional range of the "Memobase Plus Basic" license and also:
	 Detection and assessment of sensor condition
	 Multipoint calibration and adjustment of pH sensors
	 Monitoring of defined limits for measurements and adjustments
	• Data for Endress+Hauser reference solutions and sample data transferred via barcode
Pharmaceutics compliance	Functional range of the "Memobase Plus Basic" license and also: Advanced user administration

Demo version

A demo version can be used for free without any obligations. Neither real sensors nor MemoLink need be connected for this purpose. There are also videos available which introduce you to the software functions. For more information, please contact your Endress+Hauser service or sales office.

Accessories

The following are the most important accessories available at the time this documentation was issued.

▶ For accessories not listed here, please contact your Service or Sales Center.

Kits

Kit CYZ71D MemoLink for Memosens (incl. USB cable)

Order No. 71163002

Kit CYZ71D USB cable

Order No. 71162980

Measuring cable

Memosens laboratory cable CYK20

- For digital sensors with Memosens technology
- Product Configurator on the product page: www.endress.com/cyk20

Memosens data cable CYK10

- For digital sensors with Memosens technology
- Product Configurator on the product page: www.endress.com/cyk10



Technical Information TI00118C

Memosens data cable CYK11

- Extension cable for digital sensors with Memosens protocol
- Product Configurator on the product page: www.endress.com/cyk11



Technical Information TI00118C

Standard solutions

High-quality buffer solutions from Endress+Hauser - CPY20

The secondary buffer solutions have been referenced to primary reference material of the PTB (German Federal Physico-technical Institute) or to standard reference material of NIST (National Institute of Standards and Technology) according to DIN 19266 by a laboratory accredited by the DAkkS (German accreditation body) according to DIN 17025.

Product Configurator on the product page: www.endress.com/cpy20

Conductivity calibration solutions CLY11

Precision solutions referenced to SRM (Standard Reference Material) by NIST for qualified calibration of conductivity measuring systems in accordance with ISO 9000

- \blacksquare CLY11-A, 74 µS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz) Order No. 50081902
- \blacksquare CLY11-B, 149.6 $\mu S/cm$ (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz) Order No. 50081903
- CLY11-C, 1.406 mS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz)
 Order No. 50081904
- CLY11-D, 12.64 mS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz)
 Order No. 50081905
- CLY11-E, 107.00 mS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz)
 Order No. 50081906



Technical Information TI00162C

COY8

Zero-point gel for oxygen and disinfection sensors

- Oxygen-free and chlorine-free gel for the verification, zero point calibration and adjustment of oxygen and disinfection measuring points
- Product Configurator on the product page: www.endress.com/coy8



Technical Information TIO1244C

Sensors

Glass electrodes

Orbisint CPS11D

- pH sensor for process technology
- Optional SIL version for connecting to SIL transmitter
- With dirt-repellent PTFE diaphragm



Technical Information TI00028C

Memosens CPS31D

- pH electrode with gel-filled reference system with ceramic diaphragm
- Product Configurator on the product page: www.endress.com/cps31d



Technical Information TI00030C

Ceraliquid CPS41D

pH electrode with ceramic junction and KCl liquid electrolyte



Technical Information TI00079C

Ceragel CPS71D

pH electrode with reference system including ion trap



Technical Information TI00245C

Memosens CPS171D

- pH electrode for bio-fermenters with digital Memosens technology
- Product Configurator on the product page: www.endress.com/cps171d



Technical Information TI01254C

Orbipore CPS91D

pH electrode with open aperture for media with high dirt load



Technical Information TI00375C

Orbipac CPF81D

- Compact pH sensor for installation or immersion operation
- In industrial water and wastewater
- Product Configurator on the product page: www.endress.com/cpf81d



Technical Information TI00191C

Enamel pH electrodes

Ceramax CPS341D

- pH electrode with pH-sensitive enamel
- Meets highest demands of measuring accuracy, pressure, temperature, sterility and durability
- Product Configurator on the product page: www.endress.com/cps341d



Technical Information TI00468C

ORP sensors

Orbisint CPS12D

ORP sensor for process technology



Technical Information TI00367C

Ceraliquid CPS42D

ORP electrode with ceramic junction and KCl liquid electrolyte



Technical Information TI00373C

Ceragel CPS72D

ORP electrode with reference system including ion trap



Technical Information TI00374C

Orbipac CPF82D

- Compact ORP sensor for installation or immersion operation in process water and wastewater
- Product Configurator on the product page: www.endress.com/cpf82d



Technical Information TI00191C

Orbipore CPS92D

ORP electrode with open aperture for media with high dirt load



Technical Information TI00435C

pH ISFET sensors

Tophit CPS441D

- Sterilizable ISFET sensor for low-conductivity media
- Liquid KCl electrolyte



Technical Information TI00352C

Tophit CPS471D

- Sterilizable and autoclavable ISFET sensor for food and pharmaceutics, process engineering
- Water treatment and biotechnology



Technical Information TI00283C

Tophit CPS491D

ISFET sensor with open aperture for media with high dirt load



Technical Information TI00377C

pH and ORP combined sensors

Memosens CPS16D

- Combined pH/ORP sensor for process technology
- With dirt-repellent PTFE diaphragm
- With Memosens technology
- Product Configurator on the product page: www.endress.com/cps16D



Technical Information TI00503C

Memosens CPS76D

- Combined pH/ORP sensor for process technology
- Hygienic and sterile applications
- With Memosens technology
- Product Configurator on the product page: www.endress.com/cps76d



Technical Information TI00506C

Memosens CPS96D

- Combined pH/ORP sensor for chemical processes
- With poison-resistant reference with ion trap
- With Memosens technology
- Product Configurator on the product page: www.endress.com/cps96d



Technical Information TI00507C

Conductivity sensors with inductive measurement of conductivity

Indumax CLS50D

- High-durability inductive conductivity sensor
- For standard and hazardous area applications
- Product Configurator on the product page: www.endress.com/cls50d



Technical Information TI00182C

Indumax H CLS54D

- Inductive conductivity sensor
- With certified, hygienic design for foodstuffs, beverages, pharmaceuticals and biotechnology
- Product Configurator on the product page: www.endress.com/cls54d



Technical Information TI00508C

Conductivity sensors with conductive measurement of conductivity

Condumax CLS15D

- Conductive conductivity sensor
- For pure water, ultrapure water and hazardous area applications
- Product Configurator on the product page: www.endress.com/CLS15d



Technical Information TI00109C

Condumax CLS16D

- Hygienic, conductive conductivity sensor
- For pure water, ultrapure water and Ex applications
- With EHEDG and 3A approval
- Product Configurator on the product page: www.endress.com/CLS16d



Technical Information TI00227C

Condumax CLS21D

- Two-electrode sensor in pluq-in head version version
- Product Configurator on the product page: www.endress.com/CLS21d



Technical Information TI00085C

Memosens CLS82D

- Four-electrode sensor
- With Memosens technology
- Product Configurator on the product page: www.endress.com/cls82d



Technical Information TI01188C

Oxygen sensors

Oxymax COS22D

- Sterilizable sensor for dissolved oxygen
- With Memosens technology or as an analog sensor
- Product Configurator on the product page: www.endress.com/cos22d



Technical Information TI00446C

Oxymax COS51D

- Amperometric sensor for dissolved oxygen
- With Memosens technology
- Product Configurator on the product page: www.endress.com/cos51d



Technical Information TI00413C

Memosens COS81D

- Sterilizable, optical sensor for dissolved oxygen
- With Memosens technology
- Product Configurator on the product page: www.endress.com/cos81d



Technical Information TI01201C





