


Spare parts Promass 200 / TB2, Prosonic Flow 200, Promag 200, Prowirl 200, Levelflex FMP5x and Micropilot FMR5x - housing
Ersatzteile Promass 200 / TB2, Prosonic Flow 200, Promag 200, Prowirl 200, Levelflex FMP5x und Micropilot FMR5x - Gehäuse



Sensor	Promass 200 / TB2
Messaufnehmer	Prosonic Flow 200
Probe, antenna	Promag 200
Sonde, Antenne	Prowirl 200
	Levelflex FMP5x
	Micropilot FMR5x