

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

Endress+Hauser Flowtec AG
Kägenstrasse 7, 4153 Reinach BL, Switzerland

the Certificate No: **Z00671** for the PROFIBUS Device:

Model Name: PROMASS 83 DP
Revision: Profile 3.0; SW/FW: All 5; HW: All 5
GSD: EH_1529.GSD; File Version 0.7
PA039742.GSD

This certificate confirms that the product has successfully passed the certification tests with the following scope:

<input checked="" type="checkbox"/>	DP-V0	MS0, Sync, Freeze, Set_Slave_Add
<input checked="" type="checkbox"/>	DP-V1	MS1, MS2, I&M
<input checked="" type="checkbox"/>	Profiles	PA Devices Profile 3.0
<input checked="" type="checkbox"/>	Physical Layer	RS485

Test Report Number: **AIS 3008 PA 01/03**
Authorized Test Laboratory: **AIS, München, Germany**

The tests were executed in accordance with the following documents:
"Test Specifications for PROFIBUS DP Slaves, Version 3.0 from November 2005" and
"Test Specifications for PA Devices, Profile 3.0, Version 4.1, November 2004".

This certificate is granted according to the document:
"Framework for testing and certification of PROFIBUS and PROFINET products".

For all products that are placed in circulation by July 18, 2014 the certificate is valid for life.




(Official in Charge)

Board of PROFIBUS Nutzerorganisation e. V.



(Jörg Freitag)



(K.-P. Lindner)